

Give Your Patients Peace of Mind with the Babycom™ Home Doppler Fetal Heartbeat Monitor



As a physician, you know that for many women, pregnancy is a time of nervousness over the health and well being of their developing baby. This is especially true if a woman has complications.

Studies show that the heart rate of a developing fetus can be affected when a pregnant woman suffers anxiety and experiences stress-induced changes in her heart rate and blood pressure. (Columbia University-led study) If you have high-risk pregnancy patients that you'd like to monitor more closely or moms-to-be who are stressed because they are experiencing complications or are overly anxious and worried, the Babycom Home Doppler Fetal Heartbeat Monitor may be what your patients need.

Listening to and monitoring her baby's fetal heartbeat at home with the Babycom Home Doppler Fetal Heartbeat Monitor may help alleviate some of your patient's fears and anxiety – and can provide you with an invaluable record of the fetal heartbeat rate.

The Babycom is an advanced, simple-to-use fetal heartbeat monitor that is ***specifically designed for personal use***. The patented Babycom Home Doppler Fetal Heartbeat Monitoring Device allows women to quickly, easily, and conveniently detect the fetal heartbeat and determine the fetal heart rate themselves. It is a state-of-the-art and ***affordable*** home Doppler fetal heartbeat monitoring system.



Babycom Home Doppler Fetal Heartbeat Monitor's high-quality audio signal offers extreme clarity. The sharp and clear digital display shows the fetal heartbeat rate. Additionally, Babycom's proprietary wide-angle ultrasound beam allows for accurate results simply and safely.

Why is Babycom Home Doppler Fetal Heartbeat Monitor Different Than Most Other Doppler Fetal Heartbeat Monitors?

One of the problems with some Doppler fetal heartbeat monitors is that those monitors are designed to be used by trained medical professionals. Additionally, many of the fetal heartbeat Doppler devices on the market have a comparatively narrow beam capability, which may make it difficult for a woman to detect the very small fetal heart within the relatively large volume of the developing fetus. Often women who use these Doppler devices find it difficult to find their baby's heartbeat and may panic – resulting in a frantic call to your office.

Conversely, Babycom was specifically designed to be used at home, and the Babycom's proprietary wide-angle ultrasound beam detects fetal heartbeats with minimal positioning. The shape and design of the Babycom allows your patient to easily find her baby's heartbeat herself.

Babycom Home Fetal Doppler Heartbeat Monitor is:

- Specially designed for personal use so it's easy for your patients to use
- Suitable for use after the twelfth week of pregnancy
- Safe -- Ultrasound wave levels lower than FDA limits
- Affordably priced
- FDA cleared
- Lightweight and portable so you can take it anywhere

The Babycom Home Fetal Heartbeat Monitoring Device is *affordably* priced for your patients – only \$149.95. (Most of the other Doppler fetal heartbeat monitors on the market sell for more than double that price.)

Babycom Home Fetal Heartbeat Monitor Device Specifics:

Ultrasound Transducer

System Continuous Doppler
Oscillator Frequency (MHz) 2.5
US Intensity, max (mW/cm²) 15
FHR Range (bpm) 50 – 210

Display LCD

Physical screen size (mm) 25x12
Battery voltage 9V
Battery type recommended 9V alkaline 6LR61/6LF22
Battery life 250 operation min
Operating temperature 0 to 40oC
Transport/Storage temperature -20 to 65 degrees C
Relative humidity 5 to 95%
Dimensions (mm) 130(W) x105(D) x 220(H)
Net Weight including batteries 400g
Current Consumption 50 mA

How Can Your Patients Get the Babycom?

If you have a patient that you feel would benefit from using the Babycom, please suggest the Babycom to your patients. Your patient can place their order online by visiting www.suncoastmedical.com/babycom.cfm or call toll free 888-365-0000 x81296. Even though the Babycom is very easy to use, we do recommend that you instruct your patient on the proper use of this fetal Doppler device.

If you have any questions about the Babycom, please call us toll free 888-365-0000 x81296 or e-mail info@suncoastmedical.com. **When you call, be sure and ask about our SPECIAL BABYCOM OFFER for Physicians.**



Medical Authorization Form for Babycom Home Fetal Monitor

Please return via Fax or mail (all information is required)

Sun Coast Medical
Attn: Babycom
23875 Bennington Dr. #101
Santa Clarita, CA 91354
Fax: (661) 259-6751
Phone: 888-365-0000 x81296

PATIENT'S INFORMATION

Patient's Name: _____

Patient's Address: _____

City _____ State ____ Zip: _____

Patient's Telephone Number (if known) _____

PHYSICIAN'S INFORMATION

Authorizing Physician's Name: _____

Physician's Address: _____

City _____ State ____ Zip: _____

Physician's Phone Number: _____

Physician's Signature

Date



www.suncoastmedical.com